

MATERIALS		REQUIRED
1 sheet	$\frac{3}{16}$ x 2" x 36" med. balso.	4 strips $\frac{1}{8}$ x $\frac{7}{32}$ x 36" spruce.
2 "	$\frac{1}{8}$ x 3" x 36" "	1 " $\frac{3}{32}$ x $\frac{3}{32}$ x 36" "
1 "	$\frac{1}{4}$ x 3" x 36" hard balso	1 " $\frac{1}{8}$ x $\frac{1}{16}$ x 6" bamboo.
1 "	$\frac{3}{32}$ x 2" x 36" med. balso	1 length of 6" x $\frac{1}{8}$ dowel rod.
1 "	$\frac{1}{4}$ x 2" x 36" "	1 sheet $\frac{1}{8}$ x 3" x 3" plywood.
1 "	$\frac{1}{8}$ x 3" x 36" soft balso.	1 length 18 s.w.g. plain wire.
1 "	$\frac{1}{4}$ x 2" x 3" spruce.	1 strip 2 $\frac{1}{4}$ x $\frac{1}{4}$ x 16 s.w.g. alum.
1 "	$\frac{1}{4}$ x 2" x 3" "	$\frac{1}{2}$ of 16 s.w.g. brass tube
1 strip	$\frac{1}{4}$ x $\frac{1}{16}$ x 36" med. balso.	2 sheets lightweight tissue
1 "	$\frac{1}{4}$ x $\frac{1}{8}$ x 36" "	Lead shot for nose weight

